


# Canonical Correlation Analysis for SSVEP Detection

## STUDENTS

Hasnaa Hossam

Mohamed Hossam

Thursday, 12-9-2024



I	Objective & Overview
II	Introduction
III	System Architecture
IV	Data Preprocessing
V	Canonical Correlation Analysis (CCA)
VI	Results
VII	Future Improvements
VIII	Conclusion



- **Objective:**

To apply Canonical Correlation Analysis (CCA) for classifying EEG data in real-time using an SSVEP-based system.


- **Key Features of the Application:**

- Real-time classification of EEG data from multiple subjects.
- Uses CCA for frequency detection.
- Supports future improvements with advanced techniques.




## Overview of SSVEPs:


Steady-State Visual Evoked Potentials (SSVEPs) are continuous brain responses generated when an individual visually focuses on a stimulus flickering at a specific frequency.



## Application Flow:

1. Data Acquisition: Live EEG data input.
2. Preprocessing: Filtering and epoch extraction as described.
3. CCA-based Detection: Correlates processed EEG data with reference signals.
4. Real-time Classification: Outputs frequency class in real-time for the targeted stimulus.





## Preprocessing Steps:

### 1. Filtering:

- Band-pass filter applied to the EEG data (6Hz to 80Hz) to remove noise and non-relevant frequencies.

### 2. Artifact Removal:

- Eye movement and blink artifacts reduced by instructing subjects to avoid blinking during stimulus presentation.

### 3. Epoch Extraction:

- Time windows extracted from the EEG data around the stimulus onset with a 135ms latency correction for visual processing delays.



## Canonical Correlation Analysis (CCA)

### Standard CCA:

- Correlates multi-channel EEG with predefined sinusoidal reference signals.
- Used to identify target stimulus frequencies.

### Benefits:

- No calibration needed, robust in detecting SSVEPs.

## Why CCA:

CCA provides efficient correlation of EEG signals with sinusoidal references, making it a powerful tool for SSVEP frequency detection in brain-computer interface (BCI) applications.

## Key APP Features:

- Real-time EEG data analysis.
- CCA used for frequency identification between 9.25Hz and 14.75Hz.
- High processing efficiency, suitable for multi-user environments.







## Results & Accuracy

### Accuracy per Subject:

- Subject 4: 76.38% Accuracy
- Subject 6: 79.44% Accuracy
- Subject 7: 64.86% Accuracy
- Subject 9: 61.38% Accuracy
- Average accuracy across subjects: ~60%

## Future Improvements:

- **Use of IT-CCA:**

We plan to incorporate the IT-CCA technique, which uses individualized templates for better characterization of SSVEP signals, further improving classification accuracy.

- **Additional Channels and Harmonics:**

Exploring increased electrode coverage and more harmonics to refine detection capabilities.



**Summary:**

The APP effectively applies CCA for real-time SSVEP detection, and future improvements with IT-CCA will further enhance accuracy and efficiency.



The End

**Thank you  
for listening**

